
Deferitrin

产品编号: D57992

CAS: 239101-33-8

分子式: C₁₁H₁₁NO₄S

纯度: 99.78%

InChi Key: InChIKey=OEUFNKLKCFNLN-LLVKDONJSA-N

Smiles: O=C([C@]1(C)N=C(C2=CC=C(O)C=C2O)SC1)O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 55 mg/ml (217.16 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Deferitrin is an iron chelator. It potentially for the treatment of thalassemia and iron overload.