
Lenaldekar

产品编号: D58002

CAS: 418800-15-4

分子式: C₁₈H₁₄N₄

纯度: 99.25%

InChi Key: InChIKey=QLMFOQYHZJSYJF-CIAFOILYSA-N

Smiles: N(\N=C\C1=CNC2=C1C=CC=C2)C1=CC=CC2=CC=CN=C12

外观: 固体粉末

作用通路: Cytoskeletal Signaling; PI3K/Akt/mTOR signaling; MAPK; Tyrosine Kinase/Adaptors

作用靶点: Akt; IGF-1R; S6 Kinase

溶解性: DMSO: 55 mg/mL (192.09 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Lenaldekar inhibits T-cell expansion and autoimmune encephalomyelitis. Lenaldekar causes dephosphorylation of members of the PI3 kinase/AKT/mTOR pathway and delays sensitive cells in late mitosis.