
BX517

产品编号: D58014

CAS: 850717-64-5

分子式: C₁₅H₁₄N₄O₂

纯度: 98%

InChi Key: InChIKey=DFURSNCTQGJRRX-JYRVWZFOSA-N

Smiles: O=C(Nc1cc2c(NC(=O)/C/2=C(\c2ccc[nH]2)/C)cc1)N

外观: 固体粉末

作用通路: Cytoskeletal Signaling; PI3K/Akt/mTOR signaling

作用靶点: Akt; PDK

溶解性: DMSO: ≥27 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BX-517 is a potent and selective inhibitor of PDK1.