
CP 96021

产品编号: D58028

CAS: 139401-43-7

分子式: C₂₉H₂₁FN₄OS

纯度: 99.11%

InChi Key: InChiKey=DXTYGLAQRALLEG-LFYBBSHMSA-N

Smiles: CC=1N(C=2C(N1)=CN=CC2)C3=CC=C(OCC4=CC(/C=C/C5=NC=6C(S5)=CC=C(F)C6)=CC=C4)C=C3

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 4.93 mg/mL (10 mM), Need ultrasonic

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CP 96021 is a group of 5-HT₁ receptors and 5-HT₂ receptors antagonists that have anti-inflammatory properties and are non-self anti-ulcer compounds for the treatment of gastrointestinal tract and respiratory asthma.