
Stearyldiethanolamine

产品编号: D58037

CAS: 10213-78-2

分子式: C₂₂H₄₇NO₂

纯度: 98.16%

InChi Key: InChIKey=NDLNTMNRNCENRZ-UHFFFAOYSA-N

Smiles: OCCN(CCO)CCCCCCCCCCCCCCCCCC

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 35 mg/mL (97.9 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Stearyldiethanolamine can be used in studies about antibacterial nonwoven fabric and antibacterial freshness-keeping film.