
O4I2

产品编号: D58049

CAS: 165682-93-9

分子式: C₁₂H₁₁ClN₂O₂S

纯度: 99.25%

InChi Key: InChIKey=ULUBAPWNHROTEU-UHFFFAOYSA-N

Smiles: CCOC(=O)c1csc(n1)Nc1ccc(cc1)Cl

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: OCT

溶解性: DMSO:28.3 mg/mL (100 mM);Ethanol:14.1 mg/mL (50 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: O4I2 is an effective Oct3/4 inducer in various human cell lines including human fibroblasts.