
Benproperine phosphate

产品编号: D58053

CAS: 19428-14-9

分子式: C₂₁H₃₀NO₅P

纯度: 100%

InChi Key: InChiKey=MCVUURBOSHQXMK-UHFFFAOYSA-N

Smiles: OP(O)(O)=O.CC(COc1ccccc1Cc1ccccc1)N1CCCCC1

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: Microtubule Associated

溶解性: DMSO:60mg/mL;Water:81mg/mL;Alcohol:3mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Benproperine phosphate is a cough suppressant.