
CD73-IN-3

产品编号: D58058

CAS: 2375815-63-5

分子式: C₁₅H₁₈N₄O₂

纯度: 99.59%

InChi Key: InChIKey=WHRUIISQCORGKK-KOLCDFICSA-N

Smiles: CC1=NN=C(C=C1[C@H]2C[C@@H]2C(C)C)C3=CNC(=O)NC3=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: CD73

溶解性: DMSO: 22.5 mg/mL (78.58 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LY-3475070 is a potent and selective CD73 inhibitor with an IC₅₀ of 28 nM, and its use alone or in combination with pembrolizumab is currently being tested in patients with advanced cancer.