
LMT-28

产品编号: D58068

CAS: 1239600-18-0

分子式: C₁₇H₂₉NO₄

纯度: 99.88%

InChi Key: InChIKey=UDXWSYOXIRPYFK-ZHDDOTHNSA-N

Smiles: C([C@H])([C@@H](C(CCCCC)=C)O)C(=O)N1[C@@H]([C@H](C)C)COC1=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: IL Receptor

溶解性: DMSO : 90.0 mg/mL (289.0 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: LMT-28 is an inhibitor of IL-6 and selectively inhibits IL-6-induced phosphorylation of gp130, STAT3, and JAK2.