
HM-30181 mesylate monohydrate

产品编号: D58069

CAS: 2097125-58-9

分子式: C₃₉H₄₀N₆O₁₀S

纯度: 98.81%

InChi Key: InChIKey=DGJLRUWQERZYG-UHFFFAOYSA-N

Smiles: O.COC1cc2CCN(CCc3ccc(cc3)n4nnc(n4)c5cc(OC)c(OC)cc5NC(=O)C6=CC(=O)c7ccccc7O6)Cc2cc1OC.CS(=O)(=O)O

外观: 固体粉末

作用通路: Membrane transporter/Ion channel; Neuroscience

作用靶点: P-gp

溶解性: DMSO: 55 mg/ml (70.08 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: HM-30181 mesylate monohydrate is an oral P-glycoprotein (P-gp) inhibitor, used to enhance the oral bioavailability of P-gp substrate drugs.