
Oxyfedrine L-form HCl

产品编号: D58070

CAS: 16777-42-7

分子式: C₁₉H₂₄ClNO₃

纯度: 98.66%

InChi Key: InChiKey=JMUPNNYLJGSMPK-JPJJPTBZSA-N

Smiles: C(CCN[C@H]([C@H](O)C1=CC=CC=C1)C(=O)C2=CC(OC)=CC=C2.Cl

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Adrenergic Receptor

溶解性: DMSO: 45 mg/mL (128.63 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Oxyfedrine L-form HCl is a partial agonist at beta-adrenergic receptors and acts as a coronary vasodilator and cardiotonic agent.