
PU23

产品编号: D58072

CAS: 817635-93-1

分子式: C₂₁H₁₉N₃O₃S₂

纯度: 98%

InChi Key: InChIKey=RYIRRFVWFHETH-UHFFFAOYSA-N

Smiles: O=C(c1ccccc1)NC(Nc(cc1)ccc1S(NCc1ccccc1)(=O)=O)=S

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: MRP

溶解性: DMSO: 55 mg/mL (129.25 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PU 23 is an inhibitor of multidrug-resistant protein 4 (MRP4).