
Tacrine

产品编号: D58078

CAS: 321-64-2

分子式: C₁₃H₁₄N₂

纯度: 99.94%

InChi Key: InChIKey=YLJREFDVOIBQDA-UHFFFAOYSA-N

Smiles: NC1=C(CCCC2)C2=NC3=CC=CC=C31

外观: 固体粉末

作用通路: Neuroscience

作用靶点: AChR

溶解性: DMSO: 55 mg/ml (277.41 mM); Water: Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tacrine is an indirect cholinergic agonist and centrally acting anticholinesterase. It is approved for the treatment of Alzheimer's disease.