
ML337

产品编号: D58096

CAS: 1443118-44-2

分子式: C₂₁H₂₀FNO₃

纯度: 99.09%

InChi Key: InChIKey=QBCRLDPMQHPGIM-QGZVFWFLSA-N

Smiles: C(=O)(C1=C(F)C=C(C#CC2=CC=C(OC)C=C2)C=C1)N3C[C@H](O)CCC3

外观: 固体粉末

作用通路: Neuroscience

作用靶点: GluR

溶解性: DMSO: 31.8mg/mL (90mM;)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ML337 is a negative allosteric modulator of mGlu3 with an IC₅₀ of 593 nM.