
Phosphonoformic acid trisodium salt hexa

产品编号: D58102

CAS: 34156-56-4

分子式: CH₃O₅P

纯度: 98%

InChi Key: InChIKey=ZJAOAACCNHFJAH-UHFFFAOYSA-N

Smiles: O=C(P(O)(O)=O)O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO<1.26 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Phosphonoformic acid trisodium salt hexa is an antiviral drug for the treatment of CMV retinitis.