
Methenamine Hippurate

产品编号: D58109

CAS: 5714-73-8

分子式: C₁₅H₂₁N₅O₃

纯度: 96.75%

InChi Key: InChIKey=ROAIXOJGRFKICW-UHFFFAOYSA-N

Smiles: C1N2CN3CN1CN(C2)C3.OC(=O)CNC(=O)c1ccccc1

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: Antibacterial

溶解性: DMSO:14 mg/mL (43.83 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Methenamine Hippurate is a component of Hiprex, with antibacterial activity.