
Sildenafil Mesylate

产品编号: D58115

CAS: 1308285-21-3

分子式: $C_{22}H_{30}N_6O_4S \cdot xCH_4O_3S$

纯度: 100%

InChi Key: InChIKey=WEWNUXJEVSROFW-UHFFFAOYSA-N

Smiles: CS(O)(=O)=O.CCCC1=NN(C)C2=C1N=C(NC2=O)C1=CC(=CC=C1OCC)S(=O)(=O)N1CCN(C)CC1

外观: 固体粉末

作用通路: Metabolism

作用靶点: PDE

溶解性: DMSO:93 mg/mL; Ethanol:4 mg/mL; Water:92 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Sildenafil Mesylate is a mesylate form of Sildenafil, an inhibitor of Phosphodiesterase 5.