
σ1 Receptor antagonist-1

产品编号: D58123

CAS: 1204401-49-9

分子式: C₂₁H₂₃NOS

纯度: 99.68%

InChi Key: InChIKey=XWMHOC DNGHEUMF-UHFFFAOYSA-N

Smiles: O=C1CC(C2CCN(CC3=CC=CC=C3)CC2)SC4=CC=CC=C14

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: Sigma receptor

溶解性: DMSO: 90 mg/mL (266.68 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: σ1 Receptor antagonist-1 is a selective sigma 1 receptor antagonist.