
GSK 2830371

产品编号: D58138

CAS: 1404456-53-6

分子式: C₂₃H₂₉ClN₄O₂S

纯度: 99.3%

InChi Key: InChIKey=IVDUVEGCMXCMSO-FQEVSTJZSA-N

Smiles: Cc1ncc(Cl)cc1NCc1ccc(s1)C(=O)N[C@@H](CC1CCCC1)C(=O)NC1CC1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Phosphatase

溶解性: DMSO:46.1 mg/mL (100 mM);Ethanol:23.1 mg/mL (50 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GSK2830371 is an orally active, allosteric Wip1 phosphatase inhibitor with IC₅₀ of 6 nM.