
(E)-Akt inhibitor-IV

产品编号: D58139

CAS: 959841-49-7

分子式: C₃₁H₂₉N₄S

纯度: 99.01%

InChi Key: InChIKey=FIKDMAZMOIQHOO-ANVLNOONSA-N

Smiles: CCN(C1/C=C/N(C)C2CCCCC2)C(CC(CC2)-C3NC(CCCC4)C4S3)C2N1C1CCCCC1.I

外观: 固体粉末

作用通路: Cytoskeletal Signaling; PI3K/Akt/mTOR signaling

作用靶点: Akt

溶解性: DMSO : 50 mg/mL (81.36 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: (E)-Akt inhibitor-IV ((E)-AKTIV) is an inhibitor of PI3K-Akt.