

---

## Pomisartan

产品编号: D58145

CAS: 144702-17-0

分子式: C<sub>31</sub>H<sub>30</sub>N<sub>4</sub>O<sub>2</sub>

纯度: 99.51%

InChi Key: InChIKey=SOYCBUUKISWFER-UHFFFAOYSA-N

Smiles: O=C(O)C=1C=CC=CC1C=2C=CC(=CC2)CN3C4=CC(=CC(=C4N=C3CC)C)C=5N=C6N(C5)CCCC6

外观: 固体粉末

作用通路: Endocrinology/Hormones

作用靶点: RAAS

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Pomisartan is a small molecule angiotensin type 1 receptor (AT1R) antagonist that is used in the treatment of cardiovascular disease, studying heart failure and hypertension