

---

## 5-HT7 agonist 1

产品编号: D58148

CAS: 334974-31-1

分子式: C<sub>19</sub>H<sub>20</sub>ClN<sub>3</sub>

纯度: 98.61%

InChi Key: InChIKey=IKHUKLQDWXAECD-UHFFFAOYSA-N

Smiles: ClC=1C=CC=CC1CN2CCN(C=3C=CC=C4NC=CC43)CC2

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 45 mg/mL (138.1 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 5-HT7 agonist 1 is a selective agonist of 5-HT7 (IC<sub>50</sub> = 222.93 nM). 5-HT7 agonist 1 can be used in studies about CNS disorders.