
K-975

产品编号: D58155

CAS: 2563855-03-6

分子式: C₁₆H₁₄ClNO₂

纯度: 98%

InChi Key: InChIKey=KPXAHQSIROUBPH-UHFFFAOYSA-N

Smiles: C=CC(NC1=CC=C(C)C(OC2=CC=C(Cl)C=C2)=C1)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO : 215 mg/mL (747.20 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: K-975 is a highly selective, orally active TEAD inhibitor, effectively inhibiting protein-protein interactions between TEAD and YAP1/TAZ.