

---

## NG-Nitroarginine methyl ester

产品编号: D58169

CAS: 50903-99-6

分子式: C<sub>7</sub>H<sub>15</sub>N<sub>5</sub>O<sub>4</sub>

纯度: 99.72%

InChi Key: InChIKey=KCWZGJVSDFYRIX-YFKPBYRVSA-N

Smiles: COC([C@H](CCCNC(N[N+](O-)=O)=N)N)=O

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: NO Synthase

溶解性: DMSO: 2.33 mg/ml (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: NG-Nitroarginine methyl ester is a non-selective nitric oxide synthase inhibitor. It has been utilized experimentally to induce hypertension.