
CCR3 antagonist 1

产品编号: D58170

CAS: 879399-82-3

分子式: C₁₉H₂₁Cl₂N₃O₄S₂

纯度: 99.85%

InChi Key: InChIKey=QBXNXHODSIQDRT-AWEZNPCLSA-N

Smiles: C(N1C[C@H](CNC(CSC2=NC(CC(O)=O)=CS2)=O)OCC1)C3=CC(Cl)=C(Cl)C=C3

外观: 固体粉末

作用通路: Immunology/Inflammation;Microbiology/virology

作用靶点: CCR

溶解性: DMSO: 55 mg/mL (112.15 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CCR3 antagonist 1 is a potent CCR3 antagonist, used for the research of inflammatory and immunologic diseases.