
Dideoxyadenosine

产品编号: D58175

CAS: 4097-22-7

分子式: C₁₀H₁₃N₅O₂

纯度: 98%

InChi Key: InChIKey=WVXRAFOPTSTNLL-NKWVEPMBSA-N

Smiles: Nc1ncnc2n(cnc12)[C@H]1CC[C@@H](CO)O1

外观: 固体粉末

作用通路: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 2',3'-Dideoxyadenosine is an inhibitor of HIV replication with antiretroviral activity and antiviral efficacy.