

---

## Besifovir

产品编号: D58178

CAS: 441785-25-7

分子式: C<sub>10</sub>H<sub>14</sub>N<sub>5</sub>O<sub>4</sub>P

纯度: 98%

InChi Key: InChIKey=KDNSSKPZBDNJDF-UHFFFAOYSA-N

Smiles: Nc1ncc2ncn(CC3(CC3)OCP(O)(O)=O)c2n1

外观: 固体粉末

作用通路: Microbiology/virology

作用靶点: HBV

溶解性: DMSO:250 mg/mL (835.51 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Besifovir, a parent drug converted by LB80380, further metabolizes to its active form, LB80317. LB80380 is potent antiviral agent against hepatitis B virus (HBV).