
Acoramidis hydrochloride

产品编号: D58183

CAS: 2242751-53-5

分子式: C₁₅H₁₈ClFN₂O₃

纯度: 99.57%

InChi Key: InChIKey=MGFZEARHINUOMX-UHFFFAOYSA-N

Smiles: Cl.O=C(O)C1=CC=C(F)C(OCCCC=2C(=NNC2C)C)=C1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 58.5mg/mL (177.9mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Acoramidis hydrochloride is an orally active and selective kinetic stabilizer of WT and V122I-TTR (transthyretin) . Acoramidis hydrochloride is used in the study for transthyretin amyloidosis.