
HDAC8-IN-20a

产品编号: D58206

CAS: 1884231-52-0

分子式: C₁₅H₁₅NO₄

纯度: 99.01%

InChi Key: InChIKey=IGVRDNCZCPXENL-UHFFFAOYSA-N

Smiles: O=C(NO)C1=CC=C(OC)C(OCC=2C=CC=CC2)=C1

外观: 固体粉末

作用通路: Chromatin/Epigenetic;DNA Damage/DNA Repair

作用靶点: HDAC

溶解性: DMSO: 22.5 mg/mL (82.33 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: HDAC8-IN-20a is a potent and selective HDAC8 inhibitor with an IC₅₀ of 27 nM. HDAC8-IN-20a blocks activation of growth receptor survival signals.