
APNEA

产品编号: D58211

CAS: 89705-21-5

分子式: C₁₈H₂₂N₆O₄

纯度: 99.58%

InChi Key: InChiKey=XTPOZVLRZZIEBW-SCFUHWHPA-N

Smiles: N(CCC1=CC=C(N)C=C1)C2=C3C(N(C=N3)[C@@]4(O[C@H](CO)[C@@H](O)[C@H]4O)[H])=NC=N2

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Adenosine Receptor

溶解性: DMSO: 200 mg/mL (517.58 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: APNEA is a non-selective agonist of adenosine A₃ receptor.