
Betahistine mesylate

产品编号: D58214

CAS: 54856-23-4

分子式: $C_8H_{12}N_2 \cdot 2(CH_4O_3S)$

纯度: 98%

InChi Key: InChIKey=ZBJJDYGCNTNTH-UHFFFAOYSA-N

Smiles: CS(=O)(=O)O.CNCCc1ccccc1.CS(=O)(=O)O

外观: 固体粉末

作用通路: GPCR/G Protein; Immunology/Inflammation; Neuroscience

作用靶点: Histamine Receptor

溶解性: DMSO:59 mg/mL(179.7 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Betahistine mesilate is a histamine analog and H1 receptor agonist that serves as a vasodilator. It is used in MENIERE DISEASE and in vascular headaches but may exacerbate bronchial asthma and peptic ulcers.