
Haspin-IN-3

产品编号: D58215

CAS: 2416569-95-2

分子式: C₁₆H₁₀N₂O₃

纯度: 98.78%

InChi Key: InChIKey=VHCZWXXMSGKCKC-NSIKDUERSA-N

Smiles: C/C=C1C=CC(NC1)=NC=CC2=C/C3\OC=4C(C3=O)=CC=C(O)C4

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: Haspin Kinase

溶解性: DMSO : 22.5 mg/mL (80.9 mM; ultrasonic and warming and heat to 60°C)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Haspin-IN-3 is a potent haspin inhibitor with an IC₅₀ of 14 nM. Haspin-IN-3 has anticancer activity.