
Saroglitazar

产品编号: D58219

CAS: 495399-09-2

分子式: C₂₅H₂₉NO₄S

纯度: 98%

InChi Key: InChIKey=MRWFZSLZNUJVQW-DEOSSOPVSA-N

Smiles: C(COC1=CC=C(C[C@H](OCC)C(O)=O)C=C1)N2C(=CC=C2C)C3=CC=C(SC)C=C3

外观: 固体粉末

作用通路: DNA Damage/DNA Repair;Metabolism

作用靶点: PPAR

溶解性: DMSO: 22.5mg/mL (51.2mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Saroglitazar is an agonist of PPAR with EC₅₀ values of 0.65 pM and 3 nM for PPAR α and PPAR γ .