

---

## N-Methylisatoic anhydride

产品编号: D58235

CAS: 10328-92-4

分子式: C<sub>9</sub>H<sub>7</sub>NO<sub>3</sub>

纯度: 95.27%

InChi Key: InChIKey=KJMRWDHBVCNLTQ-UHFFFAOYSA-N

Smiles: CN1C(=O)OC(=O)C2=CC=CC=C12

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 60 mg/mL (338.68 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: N-Methylisatoic anhydride is a 2'-OH selective acylation agent of RNAs. It is also widely used for resolving secondary RNA structures using the SHAPE technology.