
Topotecan

产品编号: D50548

CAS: 123948-87-8

保存条件: Store at -20°C for 3 year (Powder) ; In DMSO or others solvent store at 2-4°C for two weeks,
at -20°C for six months.

产品介绍: CAS号: 123948-87-8

分子式: $C_{23}H_{23}N_3O_5$

分子量: 421.45

MDL: MFCD00870670

纯度: $\geq 99\%$

外观/性状: 粉末

溶解性: Soluble in water (200 mg/mL), and DMSO (95 mg/mL). Insoluble in ethanol.

作用靶点: Topoisomerase; Autophagy;

作用通路: Cell Cycle/DNA Damage; Autophagy;

描述: 拓扑替康 (Topotecan) 是拓扑异构酶I抑制剂。