
Tozasertib

产品编号: D50554

CAS: 639089-54-6

保存条件: Store at -20°C for 3 year (Powder) ; In DMSO or others solvent store at 2-4°C for two weeks,
at -20°C for six months.

产品介绍: CAS号: 639089-54-6

分子式: $C_{23}H_{28}N_8OS$

分子量: 464.59

MDL: MFCD13185152

纯度: $\geq 99.5\%$ (TLC)

外观: 粉末

溶解性: 溶于DMSO (≥ 106.67 mg/mL)

作用靶点: Aurora Kinase; Autophagy;

作用通路: Cell Cycle/DNA Damage; Epigenetics; Autophagy;

产品描述: 陶扎色替 (Tozasertib; VX 680; MK-0457) 是 Aurora A/B/C 激酶抑制剂, Ki 值分别为 0.6, 18, 4.6 nM。