
Btk inhibitor 2

产品编号: D58268

CAS: 1558036-85-3

分子式: C₂₄H₂₅N₅O₃

纯度: 99.49%

InChi Key: InChIKey=WTLRSSATGYETCG-QGZVFWFLSA-N

Smiles: C=CC(=O)N1CCC[C@H](C1)N2C(=C(C(=N2)C3=CC=C(C=C3)OC4=CC=CC=C4)C(=O)N)N

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: BTK

溶解性: DMSO:75 mg/mL (173.82 mM; Need ultrasonic and warming)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Btk inhibitor 2 is a BTK inhibitor extracted from patent US 20170224688 A1 (BGB-3111 analog).