
Afimetoran

产品编号: D58274

CAS: 2171019-55-7

分子式: C₂₆H₃₂N₆O

纯度: 98%

InChi Key: InChIKey=SNFVHLQYHFQOEP-UHFFFAOYSA-N

Smiles: CC(C)c1c(-c2cn3ncnc3c(C)c2C)[nH]c2c1cc(C1CCN(CC(N)=O)CC1)cc2

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: TLR

溶解性: DMSO : 98 mg/mL (220.44 mM; ultrasonic and adjust pH to 2 with HCl)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Afimetoran is a toll-like receptor antagonist, which can be used in the research of inflammatory and autoimmune diseases [1] .