
5-HT3 antagonist 3

产品编号: D58294

CAS: 120635-47-4

分子式: C₂₀H₂₁N₃O

纯度: 99.69%

InChi Key: InChIKey=NCNFDKWULDWJDS-UHFFFAOYSA-N

Smiles: O=C1C=2C=3C=CC=C4C3N(C2CCC1CN5C=CN=C5C)CCC4

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 5.5 mg/mL (17.22 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 5-HT3 antagonist 3 (Compound 15b) is a high-affinity 5-HT3 receptor antagonist. It binds to 5-HT3 receptors in rat brain cortical membranes (K_i: 0.25 nM).