

---

## Pitstop 2

产品编号: D58344

CAS: 1419320-73-2

分子式: C<sub>20</sub>H<sub>13</sub>BrN<sub>2</sub>O<sub>3</sub>S<sub>2</sub>

纯度: 100%

InChi Key: InChIKey=CGDLWHGPJPVPDU-ATVHPVEESA-N

Smiles: O=S(C1=C2C=CC=CC2=CC=C1)(NC(S/C3=C\C4=CC=C(Br)C=C4)=NC3=O)=O

外观: 固体粉末

作用通路: Apoptosis

作用靶点: Apoptosis

溶解性: DMSO : 62.5 mg/mL (132.03 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Pitstop 2 is a inhibitor of clathrin, and inhibits clathrin-mediated endocytosis (CME) by associating with the terminal domain of clathrin. It has the potential for anti-cancer research.