
GY4137

产品编号: D58432

CAS: 106740-09-4

分子式: C₁₁H₁₆NO₂PS₂·C₄H₉NO

纯度: 99.81%

InChi Key: InChIKey=YZMHNNLDUWRZFW-UHFFFAOYSA-N

Smiles: SP(N1CCOCC1)(C2=CC=C(OC)C=C2)=S.N3CCOCC3

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:12.85mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GY4137 is a novel water-soluble and slow releasing H₂S donor with vasodilator and antihypertensive activity. GY4137 also exhibits anti-inflammatory and anticancer activity.