
[Leu5]-Enkephalin TFA(58822-25-6(free bas))

产品编号: D58434

CAS: 73563-78-7

分子式: C₃₀H₃₈F₃N₅O₉

纯度: 99.63%

InChi Key: InChiKey=MPCRWMAJOJIMIB-RGRVRPFLSA-N

Smiles: CC(C)C[C@H](NC(=O)[C@H](Cc1ccccc1)NC(=O)CNC(=O)CNC(=O)[C@@H](N)Cc1ccc(O)cc1)C(=O)O.O=C(O)C(F)(F)F

外观: 固体粉末

作用通路: Endocrinology/Hormones; GPCR/G Protein; Neuroscience; Metabolism

作用靶点: Endogenous Metabolite; Opioid Receptor

溶解性: DMSO: 22.5 mg/mL (33.6 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: [Leu5]-Enkephalin is the most potent δ-receptor ago