
XY101

产品编号: D58435

CAS: 2349368-16-5

分子式: C₂₅H₂₀F₇NO₄S

纯度: 98.81%

InChi Key: InChIKey=AUIAOCKUNGZHV-UHFFFAOYSA-N

Smiles: O=C(NC=1C=CC(=CC1)C2=CC=C(C=C2F)C(O)(C(F)(F)F)C(F)(F)F)CC3=CC=C(C=C3)S(=O)(=O)CC

外观: 固体粉末

作用通路: Metabolism

作用靶点: ROR

溶解性: DMSO: 225mg/mL (399.3mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: XY101 is a selective, metabolically stable and orally available agonist of ROR γ inverse(IC₅₀ of 30 nM and a K_d of 380 nM).