
BI-1935

产品编号: D58456

CAS: 940954-41-6

分子式: C₂₄H₂₁F₃N₆O₃

纯度: 97.55%

InChi Key: InChIKey=KMESAVNRPDKZCQ-UHFFFAOYSA-N

Smiles: O=C(C(C=C1)=CN(CCOCC)C1=O)NC2=NC=C(N3N=C(C4=CC=CN=C4)C=C3C(F)(F)F)C=C2

外观: 固体粉末

作用通路: Metabolism

作用靶点: Epoxide Hydrolase

溶解性: DMSO:90 mg/mL (180.56mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BI-1935 is an inhibitor of soluble epoxide hydrolase (sEH). For diseases related to cardiovascular disease.