
BT173

产品编号: D58493

CAS: 2232180-74-2

分子式: C₁₈H₁₂BrN₃O₂

纯度: 98%

InChi Key: InChIKey=IZEKVHBNUXZFE-UHFFFAOYSA-N

Smiles: COC1=CC=C2C=C(C3=NC(C4=CC=C(Br)C=C4)=NO3)C=NC2=C1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BT173 is a novel inhibitor of homeodomain interacting protein kinase 2 (HIPK2), attenuating renal fibrosis through suppression of the TGF-β1/Smad3 pathway.