
TCN 213

产品编号: D58524

CAS: 556803-08-8

分子式: C₁₈H₂₄N₄O₂S₂

纯度: 100%

InChi Key: InChIKey=XBAZPYFIYYCZBO-UHFFFAOYSA-N

Smiles: O=C(CSc1nnc(NCc2ccccc2)s1)NCC1CCCCC1

外观: 固体粉末

作用通路: Neuroscience

作用靶点: NMDAR

溶解性: DMSO: 30 mg/mL (79.67 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TCN213 is an antagonist of NMDA receptor that has a selective for NR1/NR2A over NR1/NR2B