
Capadenoson

产品编号: D58525

CAS: 544417-40-5

分子式: C₂₅H₁₈ClN₅O₂S₂

纯度: 99.62%

InChi Key: InChIKey=CITWCLNVRIKQAF-UHFFFAOYSA-N

Smiles: N#CC1=C(C2=CC=C(OCCO)C=C2)C(C#N)=C(SCC3=CSC(C4=CC=C(Cl)C=C4)=N3)N=C1N

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Adenosine Receptor

溶解性: DMSO:50 mg/mL (96.15 mM);Water:Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Capadenoson is a selective adenosine-A1 receptor agonist.