
CC-671

产品编号: D58531

CAS: 1618658-88-0

分子式: C₂₈H₂₈N₆O₄

纯度: 98.65%

InChi Key: InChIKey=CWJLAVRXVFHDSJ-UHFFFAOYSA-N

Smiles: CNC(=O)c1ccc(Nc2nc(OC3CCCC3)c3c(c[nH]c3n2)-c2ccc3nc(C)oc3c2)c(OC)c1

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint

作用靶点: CDK

溶解性: DMSO:60 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CC-671 is a dual TTK protein kinase (IC₅₀: 0.005 μM) /CLK2 (IC₅₀: 0.006 μM) inhibitor.