

---

## CINPA1

产品编号: D58550

CAS: 102636-74-8

分子式: C<sub>23</sub>H<sub>29</sub>N<sub>3</sub>O<sub>3</sub>

纯度: 99.84%

InChi Key: InChIKey=AYQBYSPEGRYKFA-UHFFFAOYSA-N

Smiles: O=C(OCC)NC1=CC=C2C(=C1)N(C(=O)CN(CC)CC)C=3C=CC=CC3CC2

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 200 mg/mL (505.7 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CINPA1 is a selective inhibitor of constitutive androstane receptor (CAR) with an IC<sub>50</sub> of 70 nM for CAR-mediated transcription. CINPA1 can be used in studies about CAR function.