

---

## ASP-4058

产品编号: D58573

CAS: 952565-91-2

分子式: C<sub>19</sub>H<sub>12</sub>F<sub>6</sub>N<sub>4</sub>O<sub>2</sub>

纯度: 99.55%

InChi Key: InChIKey=NJNXCJPSMWKXHO-VIFPVBQESA-N

Smiles: O([C@H](C(F)(F)F)C)C1=C(C(F)(F)F)C=C(C=C1)C2=NC(=NO2)C=3C=C4C(=CC3)N=CN4

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: S1P Receptor;LPL Receptor

溶解性: DMSO : 45.0 mg/mL (101.7 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ASP-4058 is a selective, safe, and orally active second-generation agonist of sphingosine phosphate receptors 1 and 5 (S1P1 and S1P5).ASP-4058 ameliorates experimental autoimmune encephalomyelitis in mice.